

Experience sound
with Sebo HD hearing
instruments!

Sebo HD[™]
The High Fidelity Hearing Instrument[™]



Bringing **HD COLOR** to Sound.[™]

Ask your hearing care professional to fit you with
the Sebo HD hearing instrument today!

Sebo HD hearing instruments are...
Simply the best.

SeboTek hearing instruments are recognized around the world for their magnificent sound quality and their patented receiver in the canal design that has revolutionized the hearing aid industry. The groundbreaking Sebo HD line of instruments incorporates a host of new features designed to provide SeboTek's best in sound quality, noise reduction and speech enhancement technology.

Sebo HD helps users enjoy the finer things in life that they may have been missing... things like conversation with friends and family, the sounds of nature, or sounds needed for safety and security. Overall, Sebo HD users feel better about themselves and regain the independence and quality of life that hearing loss can take away.




The Sebo HD difference... *Superior Sound Quality.*

SeboTek is known as an industry leader for providing wide bandwidth hearing instruments that perform more like high fidelity headphones than a hearing aid.

The ability to hear a wide spectrum of sound brings life and brilliance to sound, allowing the listener to experience the clarity of voices, subtle sounds of nature, the vibrant tones of music, and the many other sounds that color daily life. For this reason SeboHD instruments provide the widest spectrum of sound reproduction in the hearing aid industry—up to 14kHz and 111 HD processing channels.


HD Sound Quality

Number of processing channels



HD Sound Quality

Bandwidth comparison



Experience hearing... *with high definition sound.*



High definition...

noise reduction, speech enhancement, and automatic features.

HD - Auto-Focus Directional Microphone

Having a conversation in a noisy environment is perhaps the most difficult listening situation a person with hearing loss will experience—often resulting in frustration. The Auto-FOCUS feature is designed to maximize amplification in front of the user to improve conversational listening and at the same time reduce noise from behind or from the side. (Not available in HD-6)

HD - Voice Priority Noise Reduction

Sebo HD instruments use an advanced digital technology that is able to monitor the users environment and recognize the difference between noise and speech. Making millions of calculations per second the instruments will automatically reduce amplification of noise and enhance the speech signal—just one more reason that SeboTek instruments perform so well. (available on select models)

Our daily lives are filled with multiple listening situations, from a quiet conversation with a loved one to a more complex environment such as a noisy restaurant or place of worship. The Sebo HD automatic environment detection feature adjusts to subtle changes in the user's acoustic environment and automatically optimizes its advanced features to maximize listening comfort.

- **Music Enhancement**—Automatically detects music and optimizes high fidelity sound.
- **Comfort in Quiet**—Automatically reduces low level noise to maximize listening comfort.
- **Comfort in Noise**—Automatically reduces loud noise to maximize listening comfort.
- **Speech in Quiet**—Automatically enhances speech amplification in a quiet environment.
- **Speech in Noise**—Automatically enhances speech amplification and reduces noise in a loud environment.
- **Telephone**—Automatically switches to telephone mode when telephone is put to the ear.

Unmatched
Comfort.

Sebo HD instruments use a patented comfort technology that results in the instruments being so comfortable that people forget they are wearing them. SeboTek's unique patented technology creates an ultra soft and gentle fit in the ear that creates an open and natural feel for the ultimate comfort.

