



Different by Design™



SEBOTEK®
HEARING SYSTEMS



REVOLUTIONARY DESIGN & TECHNOLOGY.

Different by Design.™

Hearing: it's central to every facet of life. Whether listening to your children laugh, sharing a quiet conversation with a friend, or simply enjoying music or the sounds of nature, **it is our ability to hear and participate in the world around us that brings richness to life.**

Yet, for those with hearing loss, traditional hearing aids—while beneficial—come with limitations that can be annoying. Squealing at all the wrong times. A plugged up feeling that creates a 'head in a barrel' sensation. And amplification that can magnify background noise more than speech.

The distinctive look and design of the SeboTek PAC (Post Auricular Canal) digital hearing instrument is the first indication that it is not a typical hearing aid. In fact, **SeboTek Hearing Systems has broken the mold on traditional hearing aids, using a patented, revolutionary design and technology to deliver what matters most: simple, clear, natural sound, unmatched comfort and cosmetic appeal.**

Ask your hearing care professional about the SeboTek PAC. We think you'll agree: when it comes to hearing instruments, **you've never heard sound like this before!**




EXCEPTIONAL COSMETIC APPEAL.

Virtually Invisible.



The SeboTek PAC is designed to rest discreetly in the natural ear contours. The advantage of SeboTek's design is its small size. SeboTek separates the speaker from the sound processor which greatly reduces size and visibility. ***The tiny speaker disappears into your ear and the miniature sound processor is virtually invisible from any angle.***

SeboTek's patented design offers a discreet, virtually invisible solution for those with severe, moderate or mild hearing loss. Best of all, because the PAC doesn't require time-consuming custom earmolds, ***you can try it in the office, and wear it home the same day!***







UNMATCHED COMFORT.

Ultra-soft, Free and Natural.

The SeboTek PAC replaces the rigid plastic ear mold of traditional hearing aids. Our patented design allows the speaker to suspend freely in the ear canal, while a flexible Ultra Soft Tip creates a soft, gently sealed fit. The result is a hearing instrument that maximizes comfort and virtually eliminates feedback (whistling), without feeling plugged up or stuffy. Experience comfort without compromise! ***The PAC offers an ultra-soft, free and natural feel that's so comfortable, you'll forget you're wearing it!***





SUPERIOR SOUND QUALITY.


High Fidelity, Acoustic Excellence.

Many people who wear hearing aids complain that the sound isn't natural—that their own voice sounds like they're "talking in a barrel" and that other sounds are "artificial and tinny."

By applying advanced principles of acoustics, **the SeboTek PAC™ delivers full spectrum audio.** The PAC reaches new heights of high fidelity audio signal processing, with **advanced digital features** designed to ensure a more true replication and amplification of sound.

The PAC combines state-of-the-art 32kHz, 20-bit digital signal processing with High Fidelity Acoustic Coupling (HiFi ACTM) and an extended bandwidth to 14,000 Hz.

Results of a clinical study in which satisfied hearing aid wearers were asked to compare the sound quality of their current instruments with the SeboTek PAC.



The Sebotek PAC uses Pro-VEST™, a programmable voice enhancement software, **designed to maximize your ability to hear speech sounds and reduce distracting background noise.** The PAC also utilizes state-of-the-art directional microphone technology that helps you hear the person speaking in front of you while, at the same time, reducing the noise behind you.

The result is brilliant, natural sound that can mean more understandable speech, all in hearing instrument that's **soft, comfortable and virtually invisible.**



DIFFERENT BY DESIGN.™

You've never heard sound like this before!



Revolutionary Design & Technology.

Exceptional Cosmetic Appeal.

Unmatched Comfort.

Superior Sound Quality.

Ask Your Hearing Care Professional to Fit You With the SeboTek PAC Today!

No time-consuming custom ear molds. No waiting weeks for your hearing aid to come back from the manufacturer.

In just a few minutes, your hearing care professional can fit you with the SeboTek PAC and you can experience its superior sound quality, unmatched comfort and cosmetic appeal for yourself!

There's no obligation. ***The SeboTek PAC: Try it for yourself. We think you'll agree, when it comes to hearing instruments, you've never heard sound like this before!***

Ask your hearing care professional to fit you with
the SeboTek PAC hearing instrument today!



The benefits of hearing instruments may vary based on the type and degree of an individual's hearing loss. Consult your hearing care professional for appropriate advice.

www.SeboTek.com

SEBOTek[®]
HEARING SYSTEMS

US Patents 5,606,621; 7,139,404; 7,110,562; 7,016,512; D,501,255 S; D,529,612 S • International Patents 700,444; P19608669-6; 2,224,653; 3,811,731;
Community Design 000306733-0001 • Other Patents Pending © 2007 All Rights Reserved • BKT022-F400-A • Rev. 0711